Laboratory Supply Requisition

Facility Name	Date Requested	_ Requested by
Facility Address	ΤΘ	el #:

The supplies provided by St Joseph Hospital to your office are solely to be used for the collection, preparation, and transport of specimens which are being sent to our laboratory for testing. It is not intended that these supplies be provided for any other uses, and such a limitation is necessary to comply with applicable laws. Your acceptance of the supplies is an acknowledgement of your agreeing to the above conditions.

Item Description	Units	Quantity Requested
Affirm VPIII swab	Each	
Biohazard Specimen Bags Sm - Med - Lg (please circle size)	50/package	
Blood Culture Bottles (adult)	Set of 2	
Blood Culture (pediatric)	Each Btl	
Butterfly Collection Needle: 23G 25G (please circle size)	50/box	
Cold Pack	Each	
CT-NG Collection Kit – vaginal/endocervical, swab	Each	
CT-NG Collection Kit - urine	Each	
Culturette, Anaerobic (Grey Cap)	Each	
Culturette swabs, double (Red Cap)	50/bag	
Culturette swabs, single (White Cap)	50/bag	
MRSA Collection Swab	Each	
Needle Hub, Venipuncture	250/bag	
<i>Needle, Venipuncture, 21G (green) 22G (black)</i> (please circle size)	48/box	
Occult Blood collection kit, feces	Each	
PAP collection vial, ThinPrep [®]	25 / tray	
PAP Collection Device for ThinPrep [®] , Medscand Combo Kit	25 / bag	
PAP collection vial, SurePath [®]	25 / tray	
PAP Collection Device for SurePath [®] , Cervex-Brush (blue)	25 / bag	
Pertussis Collection Swab (must be refrigerated)	Each	
Sharps Containers - 2 Gal 4 Gal (please circle size)	Each	
Stool Culture Collection Vial (Orange)	Each	
Stool for Giardia/Crypto Collection Vials (Pink/Green)	Each Set	
Toilet Specimen Collection devices ("Hats")	Each	
Tourniquets, non-latex	250/Box	
Tube, ACD (Yellow) 5.0 mL	Each	
Tube, SST (Gold)3.5 mL4.0 mL(please circle size)	100 / flat	
Tube, Clot (Red)4.0 mL6.0 mL(please circle size)	100 / flat	
<i>Tube, Coag (Light Blue)</i> 1.8 mL 2.7 mL (please circle size)	100 / flat	
<i>Tube, EDTA (purple)</i> 2.0 mL 4.0 mL (please circle size)	100 / flat	
Tube, Gray, 2.0 mL	Each	
<i>Tube, Light Green, PST</i> 3.0 mL 4.5 mL (please circle size)	100 / flat	
Tube, Trace Metals Analysis (Dark Blue) 6.0 mL	Each	
Urine Container, 24 hour collection	Each	
Urine Specimen Preservation Kit, BD #364957 (for UA & Culture)	Each	
Urine Cups, Sterile	100/bag	
Viral Transport Media (M4/M5)	Each	
Warm packs	Each	
Wipe, Benzalkonium (towelettes for urine culture collection)	100/box	
Other: (specify)		

* Italicized Items provided only to sites that collect blood specimens

Fax completed form to 603-595-3056